

CROS/BiCROS X



CROS/BiCROS X

**HEAR WHAT  
MATTERS  
TO YOU**

[signia-hearing.com/crosx](https://signia-hearing.com/crosx)

**signia**

Life sounds brilliant.



## SOUNDS JUST MADE FOR YOU

Your unaidable hearing loss in one ear should not stop you hearing what matters to you. With the most advanced CROS/BiCROS solutions, Signia devices ensure brilliant sound so you can hear conversations and sounds coming from both sides.

With CROS X hearing solutions it is possible to hear everything and everyone in crystal-clear quality. Wherever you go.


Whether you're looking for rechargeability, Bluetooth® streaming or discreet design, Signia Xperience provides the optimal solution for you.

# HEARING WITH EASE

The game-changing technology of CROS X devices understands your personal hearing environment in greater detail than ever to deliver sound precisely tailored to you.

Thanks to the wireless transmission from the unaidable ear to the ear with better hearing, you can follow conversations effortlessly, even when someone is speaking to you on your unaidable side.


## CROS



Hearing aid receives sound via CROS on the unaidable side.

CROS wirelessly transmits sounds.

## BiCROS



Hearing aid receives sound via CROS on the unaidable side and amplifies as needed based on individual settings. A motion sensor in BiCROS Pure X devices ensures optimal hearing even while moving.

CROS wirelessly transmits sounds.



## CROS Silk X

- Nearly invisible in the ear for complete discretion
- Signia Xperience technology for a superior hearing experience
- Click Sleeves for a quick and easy instant-fit
- Remote control via Signia app and miniPocket

Compatible with Silk X



## CROS Pure 312 X

- Signia Xperience technology for clear speech understanding and a superior hearing experience
- Calls on your smartphone are automatically streamed directly to your hearing aids
- Listen to your favorite music or TV programs via Bluetooth
- Remote control via Signia app and miniPocket

Compatible with all Bluetooth Signia Xperience hearing aids\*

\* except Styletto X



## CROS Pure Charge&Go X

- Convenient lithium-ion rechargeability
- Signia Xperience technology for clear speech understanding and a superior hearing experience
- Calls on your smartphone are automatically streamed directly to your hearing aids
- Listen to your favorite music or TV programs via Bluetooth
- Remote control via Signia app and miniPocket

Compatible with all Bluetooth Signia Xperience hearing aids\*

\* except Styletto X



### Signia app

Enables remote control and direct streaming for hearing aids using Bluetooth

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners.

### **Legal Manufacturer**

Signia GmbH  
Henri-Dunant-Strasse 100  
91058 Erlangen  
Germany

Order No. 04163-99C1-7600 | Printed in Germany  
© 02.2020 Signia GmbH

Find us on:

[signia-hearing.com](https://www.signia-hearing.com)



[signiahearing](https://www.facebook.com/signiahearing)



[signiahearing](https://www.youtube.com/signiahearing)



[@signiahearing](https://www.instagram.com/signiahearing)