



# The power of conversation is for everyone



**signia**

Signia Integrated Xperience:  
a comprehensive portfolio for  
all levels of hearing loss

# RealTime Conversation Enhancement

Delivering unmatched clarity in noisy group settings — and now it's available for all levels of hearing loss

## Connectivity

Signia is the only manufacturer to offer Bluetooth® Classic, MFi, ASHA, and LE Audio connectivity



Pure Charge&Go BCT IX



Pure Charge&Go T IX



Styletto IX



Motion Charge&Go M IX



Motion Charge&Go SP IX



Motion Charge&Go P IX

NEW



Insio Charge&Go IX CIC



Insio IX CIC



Insio IX IIC



Insio IX ITE



Insio IX ITC (CA) / HS

NEW



Active Pro IX

NEW



Silk Charge&Go IX

## Power

Fully rechargeable BTE family for all levels of hearing loss — including the longest battery life in a rechargeable Super Power BTE hearing aid — powered by our transformative RealTime Conversation Enhancement technology

## Expanded Access

Signia expands access to hearing care by offering unique, discreet, and traditional hearing aids with advanced technology meeting people where they are in their hearing journey



# 86%

of wearers reported better performance with Signia IX<sup>1</sup>



# 24%

improvement in speech understanding<sup>2</sup>



RealTime Conversation Enhancement delivered an

# 80%

enhancement of the brain's reaction to speech sounds<sup>3</sup>

## Benefits of RealTime Conversation Enhancement

- A comfortable, immersive conversational experience without distractions
- Access to multiple speakers via the world's first multi-stream architecture<sup>4</sup> with speech that rises above the surrounding sound
- Keep up and contribute to the conversation

<sup>1</sup> Jensen, N., Samra, B., Best, S., Wilson, C., & Taylor, B. 2025 AudiologyOnline, Article 29273.

<sup>2</sup> Jensen NS, Wilson C, Kamkar Parsi H, Samra B, Hain J, Best S & Taylor B. 2024. Signia White Paper.

<sup>3</sup> Slugocki, C., Kuk, F., & Korhonen, P. 2026. Signia White Paper.

<sup>4</sup> Niels Søgaard Jensen, Katharina Mohnlein-Gilbert, Cecil Wilson, Nathan Berwick, Homayoun Kamkar Parsi, Sebastian Best, Barinder Samra, Brian Taylor. 2024. Signia White Paper.



## Pure Charge&Go IX Family

Stay connected and confident in every conversation. These small, rechargeable hearing aids keep you powered all day and connect easily to your phone, music, and TV with Bluetooth. Also includes a built-in telecoil and dual push buttons to allow quick adjustments on-the-go.

### Pure Charge&Go BCT IX


- Up to 36 hours of battery life with 5 hours of streaming
- Bluetooth Classic connectivity
- In-app support via Signia Assistant
- IP68-certified for sweat and dust resistance
- LED light instantly shows your hearing aids' power status



### Pure Charge&Go T IX

- Up to 34 hours of battery life with 5 hours of streaming
- Apple MFi, Android ASHA, LE Audio Bluetooth connectivity
- In-app support via Signia Assistant
- IP68-certified for sweat and dust resistance

#### Available Colors:



## Styletto IX

Sleek, stylish and slim, these rechargeable hearing aids give you all-day power and charging on the go. Stay connected with easy Bluetooth streaming for calls, music, and TV.

- Up to 17 hours of battery life with 5 hours of streaming
- Portable charging case holds up to 5 full charges
- Apple MFi, Android ASHA, LE Audio Bluetooth connectivity
- In-app support via Signia Assistant
- IP68-certified for sweat and dust resistance



#### Available Colors:





## Motion Charge&Go IX BTE Family

Designed for any level of hearing loss, Motion Charge&Go IX offers power and flexibility to fit your needs. Stay connected with Bluetooth streaming for calls, music, and TV, plus an integrated telecoil for easy access to looped sound systems.

- Up to 65 hours of battery life with 5 hours of streaming\*
- Apple MFi, Android ASHA, LE Audio Bluetooth connectivity
- In-app support via Signia Assistant
- IP68-certified for sweat and dust resistance
- LED Indicator light provides power status at a glance



\* Motion Charge&Go SP IX

## Insio IX CIC/IIC Family

Our most discreet, custom hearing aids feature Binaural OneMic Directionality 2.0 and Speech Stabilizer to reduce background noise and enhance conversational clarity.

### Insio Charge&Go IX CIC

- Rechargeable with a battery life of up to 35 hours
- Personalized control via the Signia App or miniPocket™ remote control
- IP68-certified for sweat and dust resistance




### Insio IX CIC/IIC

- Size 10 zinc-air battery powered
- Personalized control via the Signia App or miniPocket remote control
- Available in invisible-in-canal (IIC) and completely-in-canal (CIC) styles

### Available Colors:



### Available Faceplate Colors:



\*Not available with Insio IX IIC.



Insio IX ITC (CA)/HS

Insio IX ITE

## Insio IX ITC (CA)/HS/ITE

Custom, in-the-ear hearing aids with two microphones, ensuring superior speech enhancement in noisy settings.

- Size 312 zinc-air battery powered
- Optional integrated telecoil
- Personalized control via the Signia App or miniPocket remote control
- Available in in-the-canal (ITC), half shell (HS) and in-the-ear (ITE) styles



## Active Pro IX

The instant-fit earbud that's also a prescription hearing aid. Active Pro IX delivers an instant fit, hands-free calls, and seamless streaming — for people who don't accept limitations.

- Up to 34 hours of battery life with 5 hours of streaming
- Apple MFi, Android ASHA, and Bluetooth LE Audio connectivity
- In-app support via Signia Assistant
- IP68-certified for sweat and dust resistance



Active Earmold

## Available Faceplate Colors:



Black

Dark Brown

Mocha

Beige

Tan

Brown

## Available Colors:



Snow White

Rose Gold

Black

Brown

Mocha

Beige



## Silk Charge&Go IX

Instant fit, rechargeable completely-in-canal (CIC) featuring Binaural OneMic Directionality 2.0 and Speech Stabilizer to reduce background noise and enhance conversational clarity.

- Up to 28 hours of battery life
- Portable charger with up to 4 days of additional charges
- Personalized control via the Signia App or miniPocket remote control
- IP68-certified for sweat and dust resistance



Silk earmold

### Available Faceplate Colors:



Black



Mocha

## Signia App – you're in control

Let the easy-to-use Signia App empower you on your hearing journey:



**Signia Assistant** — Improve your listening experience with a 24/7 digital assistant that makes programming changes to your hearing aids – your personal hearing companion.



**Remote Control** — Adjust the volume and sound quality of your hearing aids to fit your preferences.




**Speech Boost** — Increase speech clarity and hear conversations more distinctly, even in noisy environments.



**My WellBeing** — Gain helpful insights into your social, physical, and hearing activities all in one place.



**Professional support** — Access to chat and remote programming services with your hearing care provider.



## Optional accessories



**TV Sound** — High quality streaming directly from your television to your hearing aids.



**StreamLine Mic** — Hands-free phone calls, direct streaming, and remote microphone functionality in premium stereo quality - ideal for one-on-one conversations.



**miniPocket™** — Discreet, convenient control of your hearing aids.

## CROS/BiCROS solutions




Signia IX offers CROS/BiCROS solutions in all performance levels and a variety of styles. Signia's unique **Skewed Directionality** provides improved speech understanding in noise.

### How CROS technology works

Signia CROS devices wirelessly transmit sound from the unaidable ear to the ear with better hearing enabling you to follow conversations effortlessly, even when someone is speaking on your unaidable side.

#### CROS


Hearing aid receives sound via CROS on the unaidable side.



CROS wirelessly transmits sounds to good ear.

#### BiCROS

Hearing aid receives sound via CROS on the unaidable side. Plus, a motion sensor in select BiCROS devices ensures optimal hearing even while moving.



CROS wirelessly transmits sounds to the better ear, which are amplified as needed based on individual settings.

## Tinnitus therapies

Signia hearing aids offer three different strategies to alleviate the effect of tinnitus, a condition where you hear sounds that are not actually there. Your Hearing Care Professional can help find the ideal solution for your needs.



### Static noise tinnitus therapy signals

Soft static therapy signals (Pink, White, Speech, High Tone, Brownian, or custom) effectively distract you from the tinnitus.



### Ocean wave tinnitus therapy signals

Ocean wave therapy signals mimic the sound of the sea to provide a positive, soothing, and stress-relieving listening experience.



### World's first: Tinnitus Notch Therapy

A clinically-proven tinnitus therapy approach that's been shown to be especially effective for people who experience tonal tinnitus, the most common type of tinnitus. Notch Therapy treats your tinnitus inaudibly and unobtrusively.

## About Signia

Signia is one of the leading hearing aid brands in the world. Our goal is to enhance human performance through iconic innovations that shape the prescription audio category. We refuse to accept hearing loss as a limitation. It's not about correcting a loss, it's about gaining an edge. Not just performing, but performing brilliantly. Not just hearing again, but hearing better than ever. Hearing every word, beat, and breath. With Signia, you can Be Brilliant.

Note: Not all features are available at all performance levels.

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by WSAUD A/S is under license. Other trademarks and trade names are those of their respective owners. Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play and the Google Play logo are trademarks of Google Inc.

Signia is a registered trademark of WSAUD A/S SI-24020



Find us on:



[www.signia.net/government-services](http://www.signia.net/government-services)



[facebook.com/signiahearing](https://facebook.com/signiahearing)



[youtube.com/signiahearing](https://youtube.com/signiahearing)



[instagram.com/signiahearing](https://instagram.com/signiahearing)

**Be  
Brilliant™**